

## CONFIDENTIAL

### Patient details

Pt ID:

Pt initials:

DoB:

Sample ID:

### Sample details

Sample date:

Sample rcd:

Delphic ID:

Report date:

### Clinic details

Clinic:

Hospital:

Town:

Physician:

**HIV-1 viral load\***

**CD4**

*absolute value and percentage*

**CD8**

*absolute value and percentage*

**CD4/CD8**

*ratio*

### VL technology

The HIV-1 viral load test has been carried out using fully automated PCR utilising the **Roche COBAS AmpliPrep** and **COBAS TaqMan**. The equipment and the HIV-1 assay carry the **CE-IVD mark** and are **FDA approved**.

COBAS, AMPLIPREP and TaqMan are trademarks of Roche

### CD4/CD8 technology

These cell counts have been carried out using the automated BD FACSCalibur™ Flow Cytometer instrument and BD Tritest™ reagents. The equipment and the assay are CE-IVD approved.

Dr C A Dale Smith  
Laboratory Director

\* TND = (<40 CPM) = Target Not Detected  
The Ct value for HIV-1 above the limit for the assay or no Ct value for HIV-1 RNA obtained  
This means that the lower limit of detection for the assay (40 CPM) has not been reached  
Assay range 40-1x10<sup>6</sup> CPM